



STERILIZATION CONTAINER SYSTEM

Standard Model — 1/1 Full Size

Series K1-10 / KP1-10 | Baskets | Accessories

The Standard Model 1/1 Full Size container system is engineered for high-capacity surgical instrument management. Available in Non-Perforated (K1-10) and Perforated (KP1-10) versions across five height options, paired with matching wire, perforated, and hybrid baskets.

PRODUCT REFERENCE

Standard Model

1/1 Full Size

Different Lids Colors



With Safety Lid




With Perforated Lid

Non Perforated		Perforated	
Code	Outside Dimension L x W x H Inside Dimension L x W x H	Code	Outside Dimension L x W x H Inside Dimension L x W x H
K1-10-110	595 x 270 x 95 mm 544 x 257 x 92 mm	KP1-10-110	595 x 270 x 95 mm 544 x 257 x 80 mm
K1-10-135	595 x 270 x 115 mm 544 x 257 x 113 mm	KP1-10-135	595 x 270 x 115 mm 544 x 257 x 101 mm
K1-10-150	595 x 270 x 140 mm 544 x 257 x 137 mm	KP1-10-150	595 x 270 x 140 mm 544 x 257 x 125 mm
K1-10-200	595 x 270 x 190 mm 544 x 257 x 187 mm	KP1-10-200	595 x 270 x 190 mm 544 x 257 x 175 mm
K1-10-260	595 x 270 x 250 mm 544 x 257 x 247 mm	KP1-10-260	595 x 270 x 250 mm 544 x 257 x 235 mm

Lids		Safety Lids	
Colors	Code	Colors	Code
Grey	K1-10-1	Grey	KS1-10-1
Yellow	K1-10-2	Yellow	KS1-10-2
Red	K1-10-3	Red	KS1-10-3
Blue	K1-10-4	Blue	KS1-10-4
Green	K1-10-5	Green	KS1-10-5
Black	K1-10-6	Black	KS1-10-6
Purple	K1-10-7	Purple	KS1-10-7
Bronze	K1-10-8	Bronze	KS1-10-8
Dimension	585 x 285 mm	Dimension	605 x 287 mm

Wire Basket



Code	Outside Dimension L x W x H
B1-5425-30	540 x 255 x 30 mm
B1-5425-50	540 x 255 x 50 mm
B1-5425-70	540 x 255 x 70 mm
B1-5425-100	540 x 255 x 100 mm

Mesh: 4x4 mm
Wire: ø 1 mm

Code	Outside Dimension L x W x H
B1-5425	510 x 240 x 5 mm

Perforated Basket



Code	Outside Dimension L x W x H
B2-5425-30	540 x 253 x 40 mm
B2-5425-50	540 x 253 x 60 mm
B2-5425-70	540 x 253 x 80 mm
B2-5425-100	540 x 253 x 110 mm

Mesh: 4x4 mm
Wire: ø 1 mm

Code	Outside Dimension L x W x H
B2-5425	540 x 256 x 5 mm

Hybrid Basket



Code	Outside Dimension L x W x H
B3-5425-50	540 x 255 x 50 mm
B3-5425-70	540 x 255 x 70 mm
B3-5425-100	540 x 255 x 100 mm

Mesh: 4x4 mm
Sides: ø 3 mm
Wire: ø 1 mm

Code	Outside Dimension L x W x H
B3-5425	540 x 255 x 5 mm

Silicone Mat



Code	Outside Dimension L x W x H
A7-5425	500 x 235 mm

Standard Model

1/1 Full Size Baskets

Wire Basket



Perforated Basket



Hybrid Basket





MARS EXPORT

MARS Global Trading Company
Global Trade | Quality | Trust

www.marsxport.com
WhatsApp: +90 546 90 90 990

STERILIZATION CONTAINERS — K1-10 / KP1-10

Rigid anodized aluminum containers in Non-Perforated and Perforated variants. Five height options accommodate the full range of surgical instrument sets.

Code	Outside Dim. L×W×H	Inside Dim. L×W×H	Code	Outside Dim. L×W×H	Inside Dim. L×W×H
NON-PERFORATED (K1-10)		PERFORATED (KP1-10)			
K1-10-110	595 x 270 x 95 mm	544 x 257 x 92 mm	KP1-10-110	595 x 270 x 95 mm	544 x 257 x 80 mm
K1-10-135	595 x 270 x 115 mm	544 x 257 x 101 mm	KP1-10-135	595 x 270 x 115 mm	544 x 257 x 101 mm
K1-10-150	595 x 270 x 140 mm	544 x 257 x 137 mm	KP1-10-150	595 x 270 x 140 mm	544 x 257 x 125 mm
K1-10-200	595 x 270 x 190 mm	544 x 257 x 187 mm	KP1-10-200	595 x 270 x 190 mm	544 x 257 x 175 mm
K1-10-260	595 x 270 x 250 mm	544 x 257 x 247 mm	KP1-10-260	595 x 270 x 250 mm	544 x 257 x 235 mm

LIDS

Grey (K1-10-1) • Yellow (K1-10-2) • Red (K1-10-3) • Blue (K1-10-4)
Green (K1-10-5) • Black (K1-10-6) • Purple (K1-10-7) • Bronze (K1-10-8)
Dimensions: 585 x 285 mm

SAFETY LIDS

Grey (KS1-10-1) • Yellow (KS1-10-2) • Red (KS1-10-3) • Blue (KS1-10-4)
Green (KS1-10-5) • Black (KS1-10-6) • Purple (KS1-10-7) • Bronze (KS1-10-8)
Dimensions: 605 x 287 mm

WIRE BASKETS — SERIES B1-5425

Mesh: 4×4 mm | Wire: ø 1 mm | Maximum steam penetration for lightweight instrument sets.

Mesh: 4×4 mm | Wire: ø 1 mm

Lid: B1-5425 | 510 x 240 x 5 mm

Code	Outside Dim. L×W×H	Lid Code	Lid Dim. L×W×H
B1-5425-30	540 x 255 x 30 mm	B1-5425	510 x 240 x 5 mm



MARS EXPORT

MARS Global Trading Company
Global Trade | Quality | Trust

www.marsxport.com
WhatsApp: +90 546 90 90 990

B1-5425-50	540 x 255 x 50 mm
B1-5425-70	540 x 255 x 70 mm
B1-5425-100	540 x 255 x 100 mm

PERFORATED BASKETS — SERIES B2-5425

Box with 10 mm feet | Mesh: 4x4 mm | Wire: ø 1 mm | Elevated design for full steam circulation.

Box with 10 mm feet | Mesh: 4x4 mm | Wire: ø 1 mm
Lid: B2-5425 | 540 x 256 x 5 mm

Code	Outside Dim. LxWxH	Lid Code	Lid Dim. LxWxH
B2-5425-30	540 x 253 x 40 mm	B2-5425	540 x 256 x 5 mm
B2-5425-50	540 x 253 x 60 mm		
B2-5425-70	540 x 253 x 80 mm		
B2-5425-100	540 x 253 x 110 mm		

HYBRID BASKETS — SERIES B3-5425

Mesh: 4x4 mm | Sides: ø 3 mm | Wire: ø 1 mm | Mesh base with solid perforated sides for enhanced retention.

Mesh: 4x4 mm | Sides: ø 3 mm | Wire: ø 1 mm
Lid: B3-5425 | 540 x 255 x 5 mm

Code	Outside Dim. LxWxH	Lid Code	Lid Dim. LxWxH
B3-5425-50	540 x 255 x 50 mm	B3-5425	540 x 255 x 5 mm
B3-5425-70	540 x 255 x 70 mm		
B3-5425-100	540 x 255 x 100 mm		



MARS EXPORT

MARS Global Trading Company
Global Trade | Quality | Trust

www.marsxport.com
WhatsApp: +90 546 90 90 990

ACCESSORIES — SILICONE MAT

SILICONE MAT

Code: AT-5425

Outside Dimension: 500 x 235 mm

Anti-slip silicone mat for secure instrument holding during sterilization, storage, and transport.

Compatible with all 1/1 Full Size containers in the AT-5425 series.

STANDARDS & COMPLIANCE

TR Origin

Türkiye

Storage

Cool & dry — away from sunlight

✓ Certifications

ISO · HACCP

GET IN TOUCH

For orders, distribution inquiries, and wholesale partnerships

www.marsxport.com | WhatsApp: +90 546 90 90 990